

PROTON IRRADIATION SITE FOR SI-DETECTORS AT THE BONN ISOCHRONOUS CYCLOTRON

D. Sauerland, R. Beck, P.D. Eversheim, HISKP, Bonn, Germany J. Dingfelder, P. Wolf, SiLab, Bonn, Germany






The Bonn Isochronous Cyclotron provides proton, deuteron, alpha particle and other light ion beams with a charge-to-mass ratio Q/A of >= 1/2 and kinetic energies ranging from 7 to 14 MeV per nucleon.

At a novel irradiation site, a 14 MeV proton beam with a diameter of a few mm is utilized to homogeneously irradiate silicon detectors, so-called devices-under-test, to perform radiation hardness studies. Homogeneous irradiation is achieved by moving the device through the beam in a row-wise scan pattern with constant velocity and a row separation smaller than the beam diameter. During the irradiation procedure, the beam parameters are continuously measured non-destructively using a calibrated, secondary electron emission-based beam monitor, installed at the exit window of the beamline. The diagnostics and the irradiation procedure ensure a homogeneous irradiation with a proton fluence error of < 2 %.


In this work, an overview of the accelerator facility is given and the irradiation site with its beam diagnostics is presented in detail.




providable ions	p, d, α ,, ¹⁶ O ⁶⁺
energy $(h = 3, Q/A \ge \frac{1}{2})$	7 to $14 \mathrm{MeV}/A$
beam current (ext.)	$\lesssim 1 \mu A$
injection / extraction radius	38 mm / 910 mm
number of revolutions	approx. 120
hill sectors	$3 \times 40^{\circ}$, 0° spiral angle
hill / valley field strength	1.9 / 0.7 T (max.)
flutter	0.62
dees	$3 \times 40^{\circ}, 40 \text{kV} \text{ (max.)}$
cyclotron harmonic h	3,9
rf frequency $v_{\rm rf}$	20.1 to 28.5 MHz
hor. / vert. emittance	16 / 22 mm mrad
relative energy width	4×10^{-3}


Position and angle stabilization of extracted beam via slit appertures

Beam Handling System: Symmetric/asymmetric double-bend monochromator or achromator



fluence distribution (row-wise)





Beam-Driven

0.3 0.4 0.5 0.6

SEMv3 sum voltage U_ Σ / V

SEM calibration

Irradiation Procedure

Pattern: Row-wise scan pattern with constant hor. velocity and a row separation smaller than beam diameter


Scan pattern is repeated until desired aim fluence is applied onto the DUT


Procedure: Autonomously and beam-driven, adapting to changing beam conditions

FC linare stages

Beam condition check (current, stability, position) at turning points, pauses until requirements are met

Results: Variance of proton fluence within the DUT region is insignificant with respect to the rel. uncertainty of the beam current measurement (<1.7%)







Extraction septum



Magnetic yoke of the cyclotron



South-West Dee

